

DAFNEBEN™

CE Class I Medical Device

8 Vaginal Capsules

Indications: Dafneben helps restore the normal endovaginal physiological condition.

The internal application of the Guar Gum gelling forms a hydrogel on the vaginal mucosa a few minutes after the dissolution of the capsules. Such hydrogel adheres to the vaginal wall and, while producing a physical barrier against pathogens, lubricates and moistens the vagina.

The probiotic Lactobacilli (Lactobacillus Crispatus, Lactobacillus Jensenii, Lactobacillus Rhamnosus, Lactobacillus Reuteri) that are scattered in the vagina rebalance the internal microflora (Doderlein complex), enhancing the barrier effect of the hydrogel and helping the restoration of the physiological pH (4.0-4.5), which is itself a protective factor.

The prebiotic inulin is a substance which is not absorbed by the body, but it represents selective nourishment for the microorganisms of the vaginal flora present in Dafneben, thus favoring their ability to colonize the vaginal ecosystem, with time reduction of restoration of the normal microflora.

The Lactobacilli are regular non-pathogenic guests of the vagina, therefore Dafneben can also be used in pregnancy and during lactation. The Guar Gum gelling agent is believed to be safe for it is widely used as a food additive (E 412).

Dafneben doesn't get absorbed, leaves no residue and is removed along with the physiological exudates normally produced by the vaginal mucosa.

Ingredients: Each hard vegetable gelatine capsule contains: Guar gum, Methylsulfonylmethane (MSM), Lactobacillus Crispatus, Lactobacillus Jensenii, Lactobacillus Rhamnosus, Lactobacillus Reuteri, Inulin, vegetable Magnesium stearate.

COMPOSITION PER CAPSULE

Lactobacillus Crispatus (10 ¹¹ ufc/g)	52,5 mg
Lactobacillus Jensenii (10 ¹¹ ufc/g)	52,5 mg
Lactobacillus Rhamnosus (10 ¹¹ ufc/g)	52,5 mg
Lactobacillus Reuteri (10 ¹¹ ufc/g)	52,5 mg
Inulin	30,0 mg

How to use:

Wash your hands thoroughly before use, extract the capsule from the housing and gently insert it as deeply as possible into the vaginal canal to facilitate the rapid dissolution, preferably in a lying position, at bedtime. The recommended dose is 1 capsule a day for 8 days, unless otherwise prescribed by your doctor.



Warnings:

Keep this product out of the sight and reach of children.

Do not take the product orally.

Do not exceed the recommended dose.

Do not use the product after the expiration date printed on the package, or if the packaging is damaged or torn.

During pregnancy and lactation use the product only after consulting your doctor.

If irritation or adverse reactions of any kind occur, consult your doctor or pharmacist.

Do not mix with other medical devices or drugs for vaginal use without first consulting your doctor.

The prolonged use of topical products applied to the surface of the body can originate phenomena of allergic sensitization. If you experience these phenomena you should stop the treatment and seek medical advice.

Contraindications:

Hypersensitivity to one or more components of the product.

Storage:

Store in a dry place away from direct light, at a temperature not exceeding 25°C.

The indicated expiration date is referred to the product as unopened and properly stored.

Do not dispose of the containers in the environment after use.

Packaging:

Cardboard pack containing a PVC / Aluminium blister with 8 capsules.

It is a CE marked medical device.



Manufactured for and Distributed by:

S&R Farmaceutici - Via dei Pioppi, 2 - 06083 Bastia Umbra (PG) - Italy
www.srfarmaceutici.com

S&R

FARMACEUTICI



Do not dispose in
the environment